

Surgical Needles

MEDITÜ

Kft. Needle Technology for Surgery.

IMPORTANT!

Dear Customers!

In as much as you would like to print the and place an order by it, please make certain that you adjust the document to 100 % when printing, you can thus avoid minor discrepancies of the potential automatic settings, which can lead to improper sizes and inappropriate order.

CATALOGUE



MEDITÚ Kft. Needle Technology for Surgery

For 30 years the Meditú Kft. has been manufacturing surgical needles and with its decades of experience it is sure to guarantee quality to its customers.

Our production facility in Makó in south-eastern Hungary offers the following range of products:

- **Surgical needles**



- **Semi-finished atraumatic needles**



- **Special, unique medical needles**

The success of our company is based on paying attention to fulfilling our customers needs and ensuring the highest quality. Our needles are manufactured from high-quality stainless steel. We place great emphasis on following all environmental compliances during the manufacturing process as well as with our products themselves.

Meditú Ltd. exports needles to over 40 countries throughout the world.

Our quality is ensured by fulfilling the requirements of the standards of ISO 9001:2015 and we circulate our products with CE conformity marking.

Diameters and Lengths of Surgical Needles (mm)

Type	Size sign	Diameter	Length
A B G AR BR GR BT GT	000	1,80	100
	00	1,60	92
	0	1,60	85
	1	1,40	80
	2	1,40	75
	3	1,30	70
	4	1,30	65
	5	1,20	60
	6	1,20	55
	7	1,10	50
	8	1,00	45
	9	1,00	40
	10	0,90	36
	11	0,80	32
	12	0,80	28
	13	0,70	25
	14	0,70	22
	15	0,60	20
	16	0,60	18
GA GAR GAT	000	2,00	90
	00	1,80	80
	0	1,80	70
	1	1,60	60
	2	1,60	55
	3	1,50	50
	4	1,40	45
	5	1,30	40
	6	1,20	35
	7	1,10	30
	8	1,00	25
E ET	000	1,30	70
	00	1,20	60
	0	1,20	55
	1	1,00	47
	2	1,00	40
	3	0,80	32
	4	0,80	28
	5	0,70	22
	6	0,60	18
OB OG OA	0	0,60	36
	1	0,60	30
	2	0,50	25
	3	0,50	21
	4	0,50	18
	5	0,50	15
H HR	00	1,40	75
	0	1,20	60
	1	1,10	52
	2	1,00	45
	3	0,90	38
	4	0,80	30
	5	0,70	22

Type	Size sign	Diameter	Length
PB PC PD PH PA	00	0,80	75
	0	0,80	65
	1	0,70	55
	2	0,70	45
	3	0,60	36
	4	0,60	30
	5	0,60	25
	6	0,60	22
	7	0,50	20
	8	0,50	18
AV	9	0,50	15
	10	0,50	12
	1	2,00	100
	2	2,00	90
	3	1,80	80
	4	1,80	70
	5	1,60	60
	6	1,60	55
	7	1,50	50
	8	1,40	45
	9	1,30	40
	10	1,20	35
	11	1,10	30
BV	12	1,00	25
	13	0,90	22
	1	2,00	100
	2	2,00	90
	3	1,80	80
	4	1,80	70
	5	1,60	60
	6	1,60	55
	7	1,50	50
	8	1,40	45
	9	1,30	40
	10	1,20	35
	11	1,10	30
BE GE BER	12	1,00	25
	13	0,90	22
	1	1,40	80
	2	1,40	75
	3	1,30	70
	4	1,30	65
	5	1,20	60
	6	1,20	55
	7	1,10	52
	8	1,10	50
	9	1,00	47
	10	1,00	45
	11	1,00	40
	12	0,90	38
	13	0,90	36
	14	0,80	34
	15	0,80	32
	16	0,80	30
	17	0,80	28
	18	0,70	25
	19	0,70	23
	20	0,70	22
	21	0,60	20
	22	0,60	18

Type	Size sign	Diameter	Length
Keith Keith-R	0	1,00	80
	1	0,90	75
	2	0,90	70
	3	0,90	65
	4	0,80	60
	5	0,80	55
	6	0,80	50
König	7	0,70	45
	1	1,20	55
	2	1,10	50
	3	1,00	45
	4	1,00	40
	5	0,90	36
	6	0,80	34
	7	0,80	32
	8	0,80	30
	9	0,80	26
	10	0,70	23
	11	0,60	20
Neumann	12	0,60	18
	1	0,70	23
Kalt	2	0,60	18
	00	0,35	25
	0	0,35	22
	1	0,35	20
	2	0,35	16
Mayo Mayo T	3	0,35	12
	4	0,35	10
	1	1,50	50
	2	1,40	45
PM/S PM/A	3	1,40	42
	4	1,30	40
	000	3,00	155
	00	3,00	150
	0	3,00	145
	1	2,00	145
	2	2,00	130
Payr G	3	2,00	115
	4	2,00	100
	5	1,80	85
	6	1,80	75
	00	0,35	25
	0	0,35	22
Payr C	1	0,35	18
	2	0,35	16
	3	0,35	12
	1	0,35	30
	2	0,35	25
	3	0,35	22
Ferguson	4	0,35	20
	5	0,35	16
	6	0,35	12
	1	1,40	74
	2	1,30	63
	3	1,10	50
	4	1,00	40
	5	0,80	34
	6	0,70	23

Catalog of surgical needles

The illustrations represent the natural size of the needles.

Eye



regular eye

– odd numbers e.g. *B 501* –



spring eye

– even numbers e.g. *B 502* –

Curvature



straight



half curved



$\frac{3}{8}$ - circle



$\frac{1}{2}$ - circle



$\frac{5}{8}$ - circle



hook-shaped



S-shaped

Point



reverse cutting



round



trocar



cutting

Order method

Send your order to the shop@meditu.hu email address.

ordering example

TYPE OF NEEDLE	SIZE SIGN
<i>B 502-</i>	<i>-11</i>













Packaging

1 doz./PE-pouch (non-sterile).
Individual packaging on request.


Special items

We produce also on request. Please contact us.









G / GR / GT $\frac{1}{2}$ - circle

Ref.-No.:	Eye	Point
G 503		
G 504		
GR 513		
GR 514		
GT 515		
GT 516		


G / GR / GT















E / ET $\frac{1}{2}$ - circle

Ref.-No.:	Eye	Point
E 529		
E 530		
ET 527		
ET 528		


E / ET















GA / GAR / GAT $\frac{1}{2}$ - circle

Ref.-No.:	Eye	Point
GA 507		
GA 508		
GAR 497		
GAR 498		
GAT 499		
GAT 500		


GA / GAR / GAT







B / BR / BT $\frac{3}{8}$ - circle

Ref.-No.:	Eye	Point
B 501		
B 502		
BR 509		
BR 510		
BT 511		
BT 512		


B / BR / BT







Keith / Keith R - straight

Ref.-No.:	Eye	Point
Keith 517		▲
Keith 518		▲
KeithR 417		●
KeithR 418		●


Keith / Keith R











H / HR $\frac{5}{8}$ - circle

Ref.-No.:	Eye	Point
H 521		▼
H 522		▼
HR 421		●
HR 422		●


H / HR







A / AR - half curved

Ref.-No.:	Eye	Point
A 505		
A 506		
AR 495		
AR 496		


A / AR







König - hook-shaped

Ref.-No.:	Eye	Point
König 525		
König 526		


König





BE / BER $\frac{3}{8}$ - circle

Ref.-No.:	Eye	Point
BE 559		▼
BE 560		▼
BER 561		●
BER 562		●


BE / BER







Ferguson $\frac{1}{2}$ - circle


Ref.-No.:	Eye	Point
Ferguson 545		●
Ferguson 546		●





Ferguson




GE / GER $\frac{1}{2}$ - circle



Ref.-No.:	Eye	Point
GE 547		▼
GE 548		▼
GER 565		●
GER 566		●

GE / GER

Mayo / Mayo T $\frac{1}{2}$ - circle


Ref.-No.:	Eye	Point
Mayo 541		▼
Mayo 542		▼
MayoT 441		▲
MayoT 442		▲

Mayo / Mayo T




PD $\frac{3}{8}$ - circle

Ref.-No.:	Eye	Point
PD 533		●
PD 534		●


PD





PB $\frac{1}{2}$ - circle

Ref.-No.:	Eye	Point
PB 537		●
PB 538		●


PB





PC - straight

Ref.-No.:	Eye	Point
PC 535		●
PC 536		●


PC





PA - half curved

Ref.-No.:	Eye	Point
PA 531		●
PA 532		●


PA





PH ⁵/₈ - circle

Ref.-No.:	Eye	Point
PH 543		●
PH 544		●


PH





OA - half curved

Ref.-No.:	Eye	Point
OA 549		▼
OA 550		▼


OA





OB ³/₈ - circle

Ref.-No.:	Eye	Point
OB 551		▼
OB 552		▼


OB




OG ¹/₂ - circle

Ref.-No.:	Eye	Point
OG 553		▼
OG 554		▼


OG





Kalt ³/₈ - circle

Ref.-No.:	Eye	Point
Kalt 557		●

Kalt





Neumann - straight

Ref.-No.:	Eye	Point
Neumann 573		●
Neumann 574		●


Neumann





Payr C - straight

Ref.-No.:	Eye	Point
Payr C 555		


Payr C





Payr G - 1/2 - circle

Ref.-No.:	Eye	Point
Payr G 655		


Payr G





PM/A - half curved

Ref.-No.:	Eye	Point
PM/A 623		


PM/A




PM/S - S- shaped

Ref.-No.:	Eye	Point
PM/S 627		


PM/S




BV - $\frac{3}{8}$ - circle

Ref.-No.:	Eye	Point
BV 580		**


BV



AV - half curved

Ref.-No.:	Eye	Point
AV 590		**

AV



Transportation occurs with the following conditions

We also accept special requests for transportation based on prior agreement, which only refers to that particular business or business partner. The validation of the contract expires upon completion of the agreement.

Prices and offers

Our prices and offers do not require engagement on our part without a written delivery agreement. For international customers, the price of delivery is always determined by the actual exchange rate at the time of the order.

Orders: occur when you give a model number, the type of the needle eye's formation and the amount.

Special manufacturing

We accept orders in small amounts as well. When submitting the offer please present a master pattern or at least one dimensional drawing or sketch. In the case of the latter two, the placer of the order should check the manufacture piece/s. As may be easily understood we are not going to start investigation if special designs are protected in some way and therefore we decline any responsibility.

Catalogue types

We reserve the right to change the construction type whenever necessary.

Delivery time

We do our best to deliver precisely, and according to the stipulated time agreed upon. Partial deliveries are also accepted by prior agreement. We reserve the right to withhold delivery or retire from the contract without granting indemnity in cases of improper business behaviour on behalf of our partners.

Dispatch

We deliver at our expense regardless of the value of the delivered products using the most favourable packaging and transport methods.

Claims for faulty goods can only be accepted in writing at the latest, two weeks after arrival of delivery.

Payment: must be effected at the term and time agreed upon based on the amount shown on the invoice. In case of delay of payment, interest rates will be applied based on actual regulations. The buyer is not permitted to withhold payments or to reckon up any service.

Property: The delivered product- even if they have been manufactured- remain our property and cannot be pledged or sold until fully paid for.

Exportation of our products requires our special agreement.

Place of completion for the delivery will be the address of the customer; according to the payment the delivering company's address and bank account.